

PRODUCT DATA SHEET

No. TCB001

1000mL Roller Bottles

Description

1000mL, Surface Treated/Non-Treated, Plug Seal Cap/Vent Cap, Sterile


Purpose

For use in applications where cell attachment is desired or not desired

Materials

Bottle: GPPS (General Polystyrene) Cap: HDPE (High-density Polyethylene) Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)

Dimensions (Unit: mm)



Features

- 2 different cap styles can be used in both open and closed systems
- Each bottle has silk-printed graduation
- Available in surface treated or non-treated
- Working volume: 100-150 mL
- Growth area: 490 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10 $^{-6}$ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



1000mL Roller Bottles

Cat. No.	Capacity	Cap Style	Surface Type	Sterilization	Qty/Pk	Qty/Case
TCB001001	1000 mL	Plug Seal Cap	Non-Treated	Yes	1	24
TCB011001	1000 mL	Plug Seal Cap	Surface Treated	Yes	1	24
TCB002001	1000 mL	Vent Cap	Non-Treated	Yes	1	24
TCB012001	1000 mL	Vent Cap	Surface Treated	Yes	1	24