

## Culture Tubes

### Description

Culture Tubes, Conical Bottom or Round Bottom, Sterile or Non-sterile

### Purpose

Culture tubes is ideal for clinical sample storage, bacterial tissue culture, and molecular biology testing consumables


### Materials

Tube: PP (Polypropylene)/PS (Polystyrene)

Cap: LDPE (Low-density Polyethylene) Color: Natural



### Dimensions (Unit: mm)



### Features

- Available with high-quality polypropylene or polystyrene
- Colorless and transparent are perfect for biological reaction experiments
- Polypropylene tubes are autoclavable at 121°C and freezable to -80°C
- Polystyrene tubes are not autoclavable
- Thicker walls will not crush
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10<sup>6</sup> (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



## Culture Tubes

Cat. No.	Capacity	Tube Material	Bottom Style	Cap Style	Sterilization	Packaging Configuration
TUB000004	4 mL	PP	Conical Bottom		No	Bulk, 1000/case
TUB010004	4 mL	PS	Conical Bottom		No	Bulk, 1000/case
TUB020004	4 mL	PP	Conical Bottom	Dual-position Cap	Yes	25/bag, 500/case
TUB012004	4 mL	PS	Conical Bottom	Dual-position Cap	Yes	25/bag, 500/case
TUB000005	5 mL	PP	Round Bottom		No	Bulk, 1000/case
TUB011005	5 mL	PS	Round Bottom		No	Bulk, 1000/case
TUB022005	5 mL	PP	Round Bottom	Plug Cap	Yes	25/bag, 500/case
TUB023005	5 mL	PS	Round Bottom	Plug Cap	Yes	25/bag, 500/case



## Culture Tubes

Cat. No.	Capacity	Tube Material	Bottom Style	Cap Style	Sterilization	Packaging Configuration
TUB025005	5 mL	PP	Round Bottom	Dual-position Cap	Yes	25/bag, 500/case
TUB028005	5 mL	PS	Round Bottom	Dual-position Cap	Yes	25/bag, 500/case
TUB000008	8 mL	PP	Round Bottom		No	Bulk, 1000/case
TUB010008	8 mL	PS	Round Bottom		No	Bulk, 1000/case
TUB020008	8 mL	PP	Round Bottom		No	125/bag, 1000/case
TUB013008	8 mL	PS	Round Bottom		No	125/bag, 1000/case
TUB002140	14 mL	PP	Round Bottom		No	Bulk, 1000/case
TUB004140	14 mL	PS	Round Bottom		No	Bulk, 1000/case
TUB012140	14 mL	PP	Round Bottom	Dual-position Cap	Yes	25/bag, 500/case
TUB014140	14 mL	PS	Round Bottom	Dual-position Cap	Yes	25/bag, 500/case