

# PRODUCT DATA SHEET



No. PPT004

# **Rack of Pipette Micro Tips**

### Description

Empty rack for 10 µL pipette micro tips, Non-sterile

### Purpose

Suitable for pipette micro tips


#### Materials


Rack: PP (Polypropylene)


### **Features**

- The high-quality polypropylene racks are fully autoclavable
- · Color-coded racks for easy volume identification
- · Hinged design for easy one-handed opening
- · Lids have a stacking feature for added stability and a secure locking tab
- · Interior support ribs prevent rack from wraping
- · Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

## Dimensions (Unit:mm)







Cat. No.	Description	Rack Color	Sterilization	Qty/Case
PPT454010	Empty rack for 10 µL pipette tips, Non-Sterile	Red	No	20







No. PPT005

# **Rack of Pipette Micro Tips**

### Description

Empty rack for 20-300 µL pipette micro tips, Non-sterile

### Purpose

Suitable for pipette micro tips


### Materials


Rack: PP (Polypropylene)


### **Features**

- The high-quality polypropylene racks are fully autoclavable
- · Color-coded racks for easy volume identification
- · Hinged design for easy one-handed opening
- · Lids have a stacking feature for added stability and a secure locking tab
- · Interior support ribs prevent rack from wraping
- · Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

## Dimensions (Unit:mm)







Cat. No.	Description	Rack Color	Sterilization	Qty/Case
PPT450200	Empty rack for 20-300 μL pipette tips, Non-Sterile	Yellow	No	20







No. PPT006

# **Rack of Pipette Micro Tips**

## Description

Empty rack for 1000 µL pipette micro tips, Non-sterile

### Purpose

Suitable for pipette micro tips


### Materials


Rack: PP (Polypropylene)


#### **Features**

- The high-quality polypropylene racks are fully autoclavable
- · Color-coded racks for easy volume identification
- · Hinged design for easy one-handed opening
- · Lids have a stacking feature for added stability and a secure locking tab
- · Interior support ribs prevent rack from wraping
- · Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

## Dimensions (Unit:mm)







Cat. No.	Description	Rack Color	Sterilization	Qty/Case
PPT450000	Empty rack for 1000 μL pipette tips, Non-Sterile	Blue	No	20









No. PPT007

# **Rack of Pipette Micro Tips**

## Description

Empty rack for 1000 µL Long pipette micro tips, Non-sterile

### Purpose

Suitable for pipette micro tips


### Materials


Rack: PP (Polypropylene)


### **Features**

- The high-quality polypropylene racks are fully autoclavable
- · Color-coded racks for easy volume identification
- · Hinged design for easy one-handed opening
- · Lids have a stacking feature for added stability and a secure locking tab
- · Interior support ribs prevent rack from wraping
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

## Dimensions (Unit:mm)







Cat. No.	Description	Rack Color	Sterilization	Qty/Case
PPT450001	Empty rack for 1000 $\mu$ L Long pipette tips, Non-Sterile	Blue	No	20

