

## Reservoir Bottles

### Description

Reservoir Bottles, PS, Sterile

### Purpose

Reservoir bottles are chemically-resistant and safe for solution storage

### Materials

Bottle: GPPS (General Polystyrene)

Cap: HDPE (High-density Polyethylene)


Color: Light green



### Features

- Available in 4 different volumes of 150, 250, 500 and 1000mL
- Light weight and heavy wall construction
- Large knurls on the reservoir bottle cap for easy screw
- Designed wide and easy access bottle mouth for efficiently and stably pour out
- Engraved graduation ensure veracity
- Special concave design for easy handling
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10<sup>-6</sup> (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

### Dimensions (Unit: mm)




1000mL



500mL



250mL



150mL

## Reservoir Bottles

Cat. No.	Capacity	Cap Type	Sterilization	Packaging Configuration
FRB000150	150 mL	GL 45	Yes	1/bag, 24/case
FRB000250	250 mL	GL 45	Yes	1/bag, 24/case
FRB000500	500 mL	GL 45	Yes	1/bag, 24/case
FRB000000	1000mL	GL 45	Yes	1/bag, 24/case