

Transfer Pipets

Description

Transfer Pipets, Clear, Sterile or Non-sterile

Purpose

Use for liquid transfer

Materials

Pipet: PE (Polyethylene)



Features

- Available with 3 different capacities of 0.2mL, 1.0mL and 3.0mL
- Transfer Pipets are perfect products for quick, safe transfer of fluid
- Long flexible stem can be bent to draw liquid from narrow or small volume tubes into the bulb
- Used in Liquid nitrogen environment
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE R

Transfer Pipets

Cat. No.	Capacity	Length	Sterilization	Packaging Configuration
PP000002	0.2 mL	68 mm	No	100/bag, 10000/case
PP000010	1.0 mL	150 mm	No	100/bag, 5000/case
PP000030	3.0 mL	155 mm	No	100/bag, 5000/case
PP003030	3.0 mL	180 mm	No	100/bag, 5000/case
PP001002	0.2 mL	68 mm	Yes	100/bag, 10000/case
PP001010	1.0 mL	150 mm	Yes	100/bag, 5000/case
PP001030	3.0 mL	155 mm	Yes	100/bag, 5000/case
PP002030	3.0 mL	180 mm	Yes	100/bag, 5000/case
PP101002	0.2 mL	68 mm	Yes	1/bag, 5000/case
PP101010	1.0 mL	150 mm	Yes	1/bag, 4000/case
PP101030	3.0 mL	155 mm	Yes	1/bag, 4000/case
PP102030	3.0 mL	180 mm	Yes	1/bag, 4000/case
PP103030	3.0 mL	180 mm	Yes	1/paper plastic bag, 200/bag, 2000/case
PP200010	1.0 mL	150 mm	No	200/bag/box, 2000/case
PP200030	3.0 mL	155 mm	No	200/bag/box, 2000/case
PP303030	3.0 mL	180 mm	No	200/bag/box, 2000/case
PP201010	1.0 mL	150 mm	Yes	1/bag, 200/box, 2000/case
PP201030	3.0 mL	155 mm	Yes	1/bag, 200/box, 2000/case
PP202030	3.0 mL	180 mm	Yes	1/bag, 200/box, 2000/case
PP205010	1.0 mL	150 mm	Yes	1/paper plastic bag, 200/box, 2000/case
PP205030	3.0 mL	155 mm	Yes	1/paper plastic bag, 200/box, 2000/case
PP203030	3.0 mL	180 mm	Yes	1/paper plastic bag, 200/box, 2000/case